

HLA-C*18
101.629-06 – including *Taq* polymerase
101.629-06u – without *Taq* polymerase
Lot No.: 0H5

Olerup SSP[®] HLA-C*18

Product number:	101.629-06 – including <i>Taq</i> polymerase 101.629-06u – without <i>Taq</i> polymerase
Lot number:	0H5
Expiry date:	2021-05-01
Number of tests:	6
Number of wells per test:	7+1

The HLA-C*18 primer set is unchanged compared to the previous *Olerup SSP[®] HLA-C*18 (Lot No. 8F3)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

C*18:01 to C*18:10, i.e. all the currently recognized HLA-C*18 alleles, will be amplified by the primers in the HLA-C*18 SSP kit; www.ebi.ac.uk/imgt/hla, 2018-July-11, release 3.33.0.

The HLA-C*18 kit enables separation of the confirmed HLA-C*18 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C*18 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

RESOLUTION IN HLA-C*18 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-C*18 SUBTYPINGS BY NON-HLA-C*18 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.