

HLA-C*18

101.629-06 – including *Taq* polymerase

101.629-06u – without *Taq* polymerase

Lot No.: 1R0

Olerup SSP[®] HLA-C*18

Product number:	101.629-06 – including <i>Taq</i> polymerase 101.629-06u – without <i>Taq</i> polymerase
Lot number:	1R0
Expiry date:	2026-04-01
Number of tests:	6
Number of wells per test:	7+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*18 LOT (9L2):

The HLA-C*18 primer set is unchanged compared to the previous *Olerup SSP[®] HLA-C*18 (Lot No. 9L2)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

C*18:01 to C*18:15, i.e. all the currently recognized HLA-C*18 alleles, will be amplified by the primers in the HLA-C*18 SSP kit; www.ebi.ac.uk/imgt/hla, 2022-January-13, release 3.47.0.

The HLA-C*18 kit enables separation of the confirmed HLA-C*18 alleles as listed in the IMGT/HLA database 3.35.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C*18 kit also enables identification of many null and alternatively expressed alleles.

RESOLUTION IN HLA-C*18 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-C*18 SUBTYPINGS BY NON-HLA-C*18 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.

