

101.611-12 – including *Taq* polymerase101.611-12u – without *Taq* polymeraseLot No.: **39F****Olerup SSP[®] HLA-Cw*03**

Product number: 101.611-12 – including *Taq* polymerase
101.611-12u – without *Taq* polymerase
Lot number: 39F
Expiry date: 2010-December-01
Number of tests: 12
Number of wells per test: 32

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*03 LOT (Y28):

Well	5'-primer	3'-primer	rationale
20	Added	Added	Primer pair added for the Cw*0347 allele.
29	Added	-	Primer added for the Cw*0350 allele.
32	New	New	New primer pair for the Cw*0348 allele.

THE NUMBER OF WELLS has been increased from 31 to 32.

ALLELE COVERAGE:

Cw*0302 to Cw*0350, i.e. all the currently recognized HLA-Cw*03 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-October-10, release 2.23.0.

RESOLUTION IN HLA-Cw*03 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*03 SUBTYPINGS BY NON-HLA-Cw*03 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.