

Olerup SSP[®] HLA-Cw*0409N

Product number:	101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase
Lot number:	19G
Expiry date:	2011-July-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*0409N Lot (74E):

The HLA-Cw*0409N primer set is unchanged compared to the previous lot.

THE PLATE LAYOUT has been changed to 2 PCR reactions in an 8-well plate.

ALLELE COVERAGE:

The Cw*0409N allele gives rise to a unique amplification pattern;
www.ebi.ac.uk/imgt/hla, 2009-May-11, release 2.25.2.

RESOLUTION IN HLA-Cw*0409N HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*0409N SUBTYPINGS BY NON-HLA-Cw*05 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.