

Olerup SSP[®] HLA-Cw*05

Product number: 101.613-12 – including *Taq* polymerase
101.613-12u – without *Taq* polymerase
Lot number: 54E
Expiry date: 2010-April-01
Number of tests: 12
Number of wells per test: 15

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*05 LOT (X61):

Well	5'-primer	3'-primer	rationale
4	-	Modified	Increased specificity of specific primer pair.
8	-	Modified	Increased specificity of specific primer pair.
11	Removed	Removed	Primer not necessary when the Cw*0502 allele has been deleted from the allele database.

THE NUMBER OF WELLS has been decreased from 16 to 15.

ALLELE COVERAGE:

Cw*0501 to Cw*0517, i.e. all the currently recognized HLA-Cw*05 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-January-11, release 2.20.0.

RESOLUTION IN HLA-Cw*05 HOMO- AND HETEROZYGOTES:

The Cw*0501,0501 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*05 SUBTYPINGS BY NON-HLA-Cw*05 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.