

HLA-Cw*14
101.625-06 – licensed for PCR
101.625-06u – not licensed for PCR
Lot No.: Y09

1

Release Note - Olerup SSP™ HLA-Cw*14 –Lot No. Y09

Product number: 101.625-06 – licensed for PCR
101.625-06u – not licensed for PCR
Lot number: Y09
Expiry date: 2009-July-01
Number of tests: 6
Number of tubes per test: 12

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*14 LOT (V93):

Tube	5'-primer	3'-primer	rationale
5	Modified	-	Increased yield of specific PCR product from the Cw*1408 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

Cw*1402 to Cw*1408, i.e. all the currently recognized HLA-Cw*14 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN HLA-Cw*14 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*14 SUBTYPINGS BY NON-HLA-Cw*14 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.