

Lot No.: **32F**

Olerup SSP[®] HLA-Cw*14

Product number: 101.625-6 – including *Taq* polymerase
101.625-6u – without *Taq* polymerase
Lot number: 32F
Expiry date: 2010-November-01
Number of tests: 6
Number of wells per test: 14

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*14 LOT (Y09):

Well	5'-primer	3'-primer	rationale
13	New	New	New primer pair for the Cw*1410 allele
14	New	New	New primer pair for the Cw*1411 allele

THE NUMBER OF WELLS has been increased from 12 to 14.

ALLELE COVERAGE:

Cw*1202 to Cw*1411, i.e. all the currently recognized HLA-Cw*12 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-October-10, release 2.23.0.

RESOLUTION IN HLA-Cw*14 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*14 SUBTYPINGS BY NON-HLA-Cw*14 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.