

## **Olerup SSP® HLA-Cw\*15**

Product number:	101.626-12 – including Taq polymerase 101.626-12u – without Taq polymerase
Lot number:	75F
Expiry date:	2011-March-01
Number of tests:	12
Number of wells per test:	16

### **CHANGES COMPARED TO THE PREVIOUS HLA-Cw\*15 LOT (20E):**

The HLA-Cw\*15 primer set is unchanged.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

Cw\*1502 to Cw\*1521, i.e. all the currently recognized HLA-Cw\*15 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-January-16, release 2.24.0.

### **RESOLUTION IN HLA-Cw\*15 HOMO- AND HETEROZYGOTES:**

The Cw\*1502,1502 genotype gives rise to a unique amplification pattern.

### **INFLUENCE ON THE INTERPRETATION OF HLA-Cw\*15 SUBTYPINGS BY NON-HLA-Cw\*15 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.