

HLA-Cw*17
101.628-06 – licensed for PCR
101.628-06u – not licensed for PCR
Lot No.: Y42

1

Release Note - Olerup SSP™ HLA-Cw*17 – Lot No. Y42

Product number: 101.628-06 – licensed for PCR
101.628-06u – not licensed for PCR
Lot number: Y42
Expiry date: 2009-September-01
Number of tests: 6
Number of tubes per test: 5

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*17 LOT (X14):

Tube	5'-primer	3'-primer	rationale
2	Modified	-	Increased specificity of specific primer pair.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

Cw*1701 to Cw*1704, i.e. all the currently recognized HLA-Cw*17 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN HLA-Cw*17 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*17 SUBTYPINGS BY NON-HLA-Cw*17 ALLELES:

No influence by non-HLA-Cw*17 alleles..

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.