

Olerup SSP[®] HLA-Cw*17

Product number: 101.628-06 – including *Taq* polymerase
101.628-06u – without *Taq* polymerase
Lot number: 69E
Expiry date: 2010-July-01
Number of tests: 6
Number of wells per test: 5

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*17 LOT (Y42):

The HLA-Cw*17 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

Cw*1701 to Cw*1704, i.e. all the currently recognized HLA-Cw*17 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-Cw*17 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*17 SUBTYPINGS BY NON-HLA-Cw*17 ALLELES:

No influence by non-HLA-Cw*17 alleles..

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.