

Olerup SSP[®] HLA-Cw*18

Product number: 101.629-06 – including *Taq* polymerase
101.629-06u – without *Taq* polymerase
Lot number: 70E
Expiry date: 2010-July-01
Number of tests: 6
Number of wells per test: 5

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*18 LOT (X86):

Well	5'-primer	3'-primer	rationale
4	Exchanged	Exchanged	New primer pair to separate the Cw*1801 and Cw*0749 alleles.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

Cw*1801 to Cw*1803, i.e. all the currently recognized HLA-Cw*18 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-April-08, release 2.21.0.

RESOLUTION IN HLA-Cw*18 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*18 SUBTYPINGS BY NON-HLA-Cw*18 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.