

## **Olerup SSP<sup>®</sup> HLA-Cw\*18**

Product number:	101.629-06 – including <i>Taq</i> polymerase 101.629-06u – without <i>Taq</i> polymerase
Lot number:	40G
Expiry date:	2011-October-01
Number of tests:	6
Number of wells per test:	5

### **CHANGES COMPARED TO THE PREVIOUS HLA-Cw\*18 LOT (70E):**

The HLA-Cw\*18 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

Cw\*1801 to Cw\*1803, i.e. all the currently recognized HLA-Cw\*18 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-October-19, release 2.27.0.

### **RESOLUTION IN HLA-Cw\*18 HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-Cw\*18 SUBTYPINGS BY NON-HLA-Cw\*18 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.