

HLA-Cw*06
101.614-12 – licensed for PCR
101.614-12u – not licensed for PCR
Lot No.: **X85**

1

Release Note - *Olerup SSP*TM HLA-Cw*06 – Lot No. **X85**

Product number: 101.614-12 – licensed for PCR
101.614-12u – not licensed for PCR
Lot number: X85
Expiry date: 2009-June-01
Number of tests: 12
Number of tubes per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-Cw*06 LOT (V92):

Tube	5'-primer	3'-primer	rationale
4	Modified	-	Increased specificity of specific primer pair.
8	-	Modified	Increased specificity of specific primer pair.
11	Modified	-	Increased specificity of specific primer pair.
14	New	New	New primer pair for the Cw*0614 allele.
15	New	New	New primer pair for the Cw*0615 allele.
16	New	New	New primer pair for the Cw*0616N allele.

THE NUMBER OF WELLS has been increased from 13 to 16..

ALLELE COVERAGE:

Cw*0601 to Cw*0616N, i.e. all the currently recognized HLA-Cw*06 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-July-09, release 2.18.0.

RESOLUTION IN HLA-Cw*06 HOMO- AND HETEROZYGOTES:

The Cw*0602,0602 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-Cw*06 SUBTYPINGS BY NON-HLA-Cw*06 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.