

KIR Genotyping
104.101-12 – licensed for PCR
104.101-12u – not licensed for PCR
Lot No.: **X70**

1

Release Note - *Olerup SSP*TM KIR Genotyping – Lot No. **X70**

Product number: 104.101-12 – licensed for PCR
104.101-12u – not licensed for PCR
Lot number: X70
Expiry date: 2009-February-01
Number of tests: 12
Number of tubes per test: 23 + 1

CHANGES COMPARED TO THE PREVIOUS KIR GENOTYPING LOT (N82):

Tube	5'-primer	3'-primer	rationale
3	-	Added	Primer added for the 2DL3*007 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

KIR alleles listed on the IPD KIR web page 2006-November-10, release 1.3.0,
www.ebi.ac.uk/ipd/kir.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.