

Olerup SSP[®] KIR Genotyping

Product number: 104.101-12 – including *Taq* polymerase
104.101-12u – without *Taq* polymerase
Lot number: 71E
Expiry date: 2010-June-01
Number of tests: 12
Number of wells per test: 23 + 1

CHANGES COMPARED TO THE PREVIOUS KIR GENOTYPING LOT (Y47):

Well	5'-primer	3'-primer	rationale
17	Exchanged	Exchanged	Primer pair exchanged to avoid co-amplification of 3DL1*009.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

KIR alleles listed on the IPD KIR web page 2008-January-18, release 2.0.0, www.ebi.ac.uk/ipd/kir.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.